

Modern Methods and Tools for Human Biosignal Analysis


Vangelis Metsis

vmetsis@txstate.edu

Computer Science – Texas State University


What is a Biosignal?

- The term *biosignal* refers to any signal that can be measured from living organisms.
- Biosignals can be classified as:
 - Bioelectrical signals (signals that originate in nerves and muscles), electrical conductance (e.g. Galvanic skin response),
 - Bioimpedance signals,
 - Bioacoustic signals,
 - Bio-optical signals (e.g. blood-oxygen saturation based on reflection or pulse rate by the change in skin color)



How are biosignals useful?

- Medical applications, e.g. heart rate, blood oxygen levels
- Sleep studies (Polysomnography)
- Affective computing, recognizing human emotion based on physiological variations.
- Human-Computer Interaction, e.g. Affect-aware Virtual Reality



Common types of Biosignals

Brain activity: Electroencephalogram (EEG) and Near-infrared spectroscopy (NIRS)

Electrical activity of muscles: Electromyogram (EMG) and Electrooculogram (EOG)

Heart activity: Heart rate (HR), and Electrocardiogram (ECG)

Electrodermal Activity (EDA), also known as skin conductance, or Galvanic Skin Response (GSR)

Air flow in and out of lungs, *Breathing effort* and rate

The *level of oxygen* in the blood (Oximetry-PPG)

How are Biosignals collected?

Invasive or non-invasive devices can be used




RS:

Photograph (PPG)

Thermometer

e

Galvanic Skin Response (GSR)



Case Study 1: Peripheral Biosignal Analysis

How are you doing?

- Humans can recognize emotional state
 - Facial expression
 - Appearance, Movement
 - Voice inflection
- Ekman and Friesen conducted research starting in the late 1960s involving facial expression as an indicator of six universal emotions: happiness, sadness, anger, fear, surprise, disgust, interest.
- Concept of affective computing described by Rosalind Picard in 1995 incorporates emotion as an input to automated systems.


A photograph of a person's hand with several sensors attached to the fingers. A large grey sensor is on the thumb, a red one on the index finger, and a black one on the middle finger. Wires connect these sensors to a black device in the bottom right corner, which has multiple colored ports (green, red, black) and a small screen. In the background, there is a laptop screen showing a colorful abstract image and some papers with handwritten notes. The title 'Beyond Human Senses' is overlaid on the left side of the image.

Beyond Human Senses

- Signals non-detectable by humans known to provide additional information regarding emotional state
- Galvanic Skin Response (Electrodermal Activity) is an example – key signal used during the administration of a lie detector test
- Technological developments have lowered the cost of monitoring physiological signals
 - *Consumer*: heart rate, movement (acceleration)
 - *Research*: brain/muscle/heart electrical activity, respiration

A bit of Biology

- Detectable emotional responses have origins in “fight or flight”
- Body attempts to stay in steady state – homeostatis
- When stressed physiological changes occur in preparation for action
 - ✓ Heart rate increases
 - ✓ Restricted peripheral blood flow
 - ✓ Sweat glands express
 - ✓ Hormones are released




Emotional Response Classification

“How does that make you feel?”


Emotional Response Classification

“How does that make you feel?”



Emotional Response Classification

“How does that make you feel?”



Experimental Design and Data Collection

Tools: BioRadio

- Physiological Data Collection
- Electrical Activity - 4/8 configurable inputs
 - Heart - electrocardiography (ECG)
 - Muscle - electromyography (EMG)
 - Eye Movement - electrooculography (EOG)
 - Brain - electroencephalography (EEG)
 - Skin Response - electrodermal activity (EDA aka GSR)
 - Breath - respiratory inductance plethysmography (RIP)
- Motion - 6-axis accelerometer, triaxial linear acceleration and angular velocity
- Auxiliary Inputs - Surface Temperature, Photoplethysmogram (PPG), Heart Rate (optical), and Oxygen Saturation (SpO₂)




Early Work VR Test Run

“Leisurely” lap

faster lap

unexpected wreck



Classification Labels – Keep it Simple!

Preliminary testing showed it was very difficult to self-classify experience

Check the box to the right that best represents how you felt during the session.

		A little							A lot						
Session	Emotion	1	2	3	4	5	6	7	1	2	3	4	5	6	7
Session 1	Excited														
	Happy														
	Sad														
	Angry														
	Scared														
	Tender														

¹ https://en.wikipedia.org/wiki/Emotion#/media/File:Emotions_-_3.png

PARTICIPANT RESPONSE

Study Title: Determination of Emotional State Through Physiological Measurement
Principal Investigator: Lee B. Hinkle
Sponsor: Dr. Vangelis Metsis

Co-Investigator/Faculty Advisor: Dr. Vangelis Metsis

1. On a scale of 1-5, how positively or negatively did you feel?

2. On a scale of 1-5, how strongly did you feel after the activity?

Activity	1	2	3	4	5
Relaxing Music	2	4			
Exciting Music	4	4			
Relaxing Video	2	4			
Exciting Video	4	5			
Tetris	4	3			
Minisoccer	4	4			

Classification Labels

Subject Self-Reporting


- Too difficult to remove headset - must be oral
- Simplified version of arousal valence model used

(1,-1) stressed	(1,0)	(1,1) excited
(0,-1)	(0,0)	(0,1)
(-1,-1) bored	(-1,0)	(-1,1) relaxed

After each session:

Question 1: "Did you find this movie exciting, relaxing, or neutral?"

Question 2: "Did you find this movie pleasant, unpleasant, or neutral?"



Sensor Setup

- Headset with headphones and Xbox game controller
- Head Radio (in back)
 - EEG/EOG/EMG(smile)
- Body Radio (in Support Belt)
 - ECG/GSR/PulseOx
 - Chest & Abdomen Respiration



Sensor List 24 total channels

Radio	Signal	Sensor Location	Qty
H_Ch1	EEG f4	High right forehead	1
H_Ch2	EOG - Horizontal	Outside of eyes	1
H_Ch3	EOG - Vertical	Above and below right eye	1
H_Ch4	EMG – Zygomaticus “smile” muscle	Right cheek	1
H_int	Accel XYZ, Gyro XYZ	Rear of head	6
B_Ch1	GSR (Electrodermal Activity)	Right index and pointer finger	1
B_Ch2	ECG	Left and right wrists	1
B_Ch3	Chest Respiration (RIP)	Chest strap	1
B_Ch4	Abdomen Respiration (RIP)	Stomach strap	1
B_Aux	Peripheral Temperature	Right pinkie finger	1
B_Aux	Heart Rate via PulseOx	Right ring finger	1
B_Aux	Blood Volume (PPG) via PulseOx	Right ring finger	1
B_Aux	Blood Oxygen (SpO2) via PulseOx	Right ring finger	1
B_int	Accel XYZ, Gyro XYZ	Right waist	6

Still Frame from Session Video

Head Radio

Oculus
Headset
Oculus
Sensor

Remote
Mic

Session View



Session Label

Head
Data

Body
Data


Quick Visual Analysis of Pendulum Swing Video Signals in BioCapture

	A	B	C	
1		MVI_004	Real Time	
2			S1_VR3H	S
3	Start H	0:17:05	3:29:41 AM	
4	Start B	0:17:21		3
5	Start of Rollercoaster [Mark]	0:17:41	3:30:15 AM	
6	Start of Drop	0:18:19	3:30:53 AM	
7	Stop Rollercoaster	0:19:31	3:32:05 AM	
8	Start Standing	0:20:03	3:32:37 AM	
9	Finish Standing	0:20:14	3:32:48 AM	
10	Jump into Unknown [Mark]	0:20:30	3:33:04 AM	
11	Pushed in VR	0:21:40	3:34:14 AM	
12	Jump into Unknown Finish	0:22:10	3:34:44 AM	
13	Ending [Mark]	0:22:16	3:34:50 AM	
14				
15	NOTE: All this done in BioCapture - need to confirm increase prec			



Classification Results

- Subject Self-Reported data for 11 segments
- Skewed toward pleasant and exciting
- No unpleasant relaxing e.g. “bored” results



Total Responses for 5 Subjects

		3
		1
		1

Introduction to Virtual Reality Demo

	1	1
	1	2

The Rose and I Movie

1	1	3

Rollercoaster 3D Movie

	1	4

Pendulum Swing 3D Movie

	1	2
	1	1

In Cell Game

1	1	
	2	1

Lost Movie

	3	
	1	1

Dream Deck Demo: Overall

	1	1
	1	1
		1

Dream Deck Demo: City Scene (Height)

		1
1	3	

Dream Deck Demo: Alien (Creepy)

2	1	1
	1	

Dream Deck Demo: Dinosaur (Scary)

		2
1	1	1


Lucky's Tale Game

4	9	17
2	11	7
0	1	4

All subjects, all sessions

Classification Results

- Subject Self-Reported data for 11 segments
- Skewed toward pleasant and exciting
- No unpleasant relaxing e.g. “bored” results



Total Responses for 5 Subjects

		3
		1
		1

Introduction to Virtual Reality Demo

	3	
	1	1

Dream Deck Demo: Overall

	1	1
	1	2

The Rose and I Movie

	1	1
	1	1
		1

Dream Deck Demo: City Scene (Height)

1	1	3

Rollercoaster 3D Movie

		1
1	3	

Dream Deck Demo: Alien (Creepy)

	1	4

Pendulum Swing 3D Movie

2	1	1
	1	

Dream Deck Demo: Dinosaur (Scary)

	1	2
	1	1

In Cell Game

		2
1	1	1

Lucky's Tale Game

1	1	
	2	1

Lost Movie

4	9	17
2	11	7
0	1	4

All subjects, all sessions

S15VR23																				
344817x29 table																				
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	
	fR	EOGH	EOGV	EMGCheek	AccelX	AccelY	AccelZ	GyroX	GyroY	GyroZ	GSR	ECG	RIPCh	RIPAb	AccelXB	AccelYB	AccelZB	GyroXB	GyroYB	
1	-0.0087	0.0234	0.0040	-0.0083	-0.0654	0.9619	-0.2239	-0.0610	-0.3662	1.0986	7.0151	0.0035	0.0014	-0.0040	0.7947	0.5400	-0.2546	1.2207	0.0610	
2	-0.0087	0.0234	0.0040	-0.0083	-0.0706	0.9573	-0.2249	0	-0.4272	1.0376	7.0151	0.0035	0.0014	-0.0040	0.7949	0.5405	-0.2454	1.0986	0.0610	
3	-0.0087	0.0234	0.0040	-0.0083	-0.0647	0.9575	-0.2231	-0.1831	-0.5493	0.9766	7.0151	0.0035	0.0014	-0.0040	0.8015	0.5322	-0.2549	1.1597	0	
4	-0.0087	0.0234	0.0040	-0.0083	-0.0674	0.9631	-0.2275	-0.1831	-0.6104	1.0376	7.0151	0.0035	0.0014	-0.0040	0.7955	0.5352	-0.2517	1.1597	-0.0610	
5	-0.0087	0.0234	0.0040	-0.0083	-0.0684	0.9578	-0.2195	-0.1221	-0.6104	1.1597	7.0151	0.0035	0.0014	-0.0040	0.7964	0.5371	-0.2517	1.0986	0	
6	-0.0087	0.0234	0.0040	-0.0083	-0.0657	0.9546	-0.2244	-0.1831	-0.4883	1.2207	7.0151	0.0035	0.0014	-0.0040	0.7930	0.5400	-0.2512	1.0986	0	
7	-0.0087	0.0234	0.0040	-0.0083	-0.0664	0.9636	-0.2236	-0.2441	-0.5493	1.2817	7.0151	0.0035	0.0014	-0.0040	0.7937	0.5374	-0.2488	1.1597	0.0610	
8	-0.0087	0.0234	0.0040	-0.0083	-0.0681	0.9590	-0.2354	-0.4272	-0.4272	1.3428	7.0124	0.0035	0.0014	-0.0040	0.7947	0.5376	-0.2410	1.2817	0.0610	
9	-0.0087	0.0234	0.0040	-0.0083	-0.0698	0.9641	-0.2295	-0.5493	-0.4272	1.6479	7.0124	0.0035	0.0014	-0.0040	0.7991	0.5317	-0.2439	1.2817	0.1831	
10	-0.0087	0.0234	0.0040	-0.0083	-0.0662	0.9629	-0.2246	-0.4883	-0.4883	1.6479	7.0124	0.0035	0.0014	-0.0040	0.7932	0.5388	-0.2451	1.2207	0.1221	
11	-0.0087	0.0234	0.0040	-0.0083	-0.0625	0.9614	-0.2246	-0.6714	-0.4272	1.7090	7.0124	0.0035	0.0014	-0.0040	0.8032	0.5371	-0.2456	1.2817	0.0610	
12	-0.0087	0.0234	0.0040	-0.0083	-0.0654	0.9578	-0.2219	-0.7324	-0.4272	1.8921	7.0124	0.0035	0.0014	-0.0040	0.8003	0.5381	-0.2463	1.2207	0.1221	
13	-0.0087	0.0234	0.0040	-0.0083	-0.0664	0.9568	-0.2224	-0.8545	-0.3662	2.0142	7.0124	0.0035	0.0014	-0.0040	0.8020	0.5388	-0.2427	1.2207	0	
14	-0.0087	0.0234	0.0040	-0.0083	-0.0632	0.9551	-0.2266	-0.9155	-0.4272	2.1973	7.0124	0.0035	0.0014	-0.0040	0.7971	0.5369	-0.2439	1.2207	-0.0610	
15	-0.0087	0.0234	0.0040	-0.0083	-0.0608	0.9587	-0.2129	-0.6714	-0.4883	2.2583	7.0124	0.0036	0.0014	-0.0040	0.7998	0.5305	-0.2446	1.2817	-0.1221	
16	-0.0087	0.0234	0.0040	-0.0083	-0.0659	0.9592	-0.2241	-0.6714	-0.4272	2.2583	7.0050	0.0036	0.0014	-0.0040	0.7983	0.5366	-0.2417	1.1597	-0.1221	
17	-0.0087	0.0234	0.0040	-0.0083	-0.0657	0.9600	-0.2231	-0.7324	-0.4272	2.4414	7.0050	0.0037	0.0014	-0.0040	0.7947	0.5393	-0.2468	1.2207	-0.2441	
18	-0.0087	0.0234	0.0040	-0.0083	-0.0630	0.9565	-0.2205	-0.6104	-0.4272	2.3804	7.0050	0.0037	0.0014	-0.0040	0.7954	0.5344	-0.2461	1.2207	-0.3052	
19	-0.0087	0.0234	0.0040	-0.0083	-0.0635	0.9578	-0.2183	-0.4883	-0.3662	2.4414	7.0050	0.0037	0.0014	-0.0040	0.7959	0.5364	-0.2512	1.2817	-0.3662	
20	-0.0087	0.0234	0.0040	-0.0083	-0.0598	0.9614	-0.2234	-0.4272	-0.4272	2.4414	7.0050	0.0037	0.0014	-0.0040	0.7959	0.5364	-0.2512	1.2817	-0.3662	
21	-0.0087	0.0234	0.0040	-0.0083	-0.0637	0.9541	-0.2283	-0.3052	-0.2441	2.5635	7.0050	0.0036	0.0014	-0.0040	0.7888	0.5352	-0.2522	1.0986	-0.0610	
22	-0.0087	0.0234	0.0040	-0.0083	-0.0559	0.9563	-0.2207	-0.3052	-0.3662	2.5024	7.0050	0.0037	0.0014	-0.0040	0.7905	0.5398	-0.2385	1.0986	-0.3052	
23	-0.0087	0.0234	0.0040	-0.0083	-0.0593	0.9575	-0.2246	-0.3662	-0.3662	2.5024	7.0050	0.0037	0.0014	-0.0040	0.7933	0.5375	-0.2419	1.2207	-0.4272	
24	-0.0087	0.0234	0.0040	-0.0083	-0.0623	0.9587	-0.2178	-0.3052	-0.4272	2.3804	7.0101	0.0037	0.0014	-0.0040	0.7883	0.5308	-0.2468	1.3428	-0.3052	
25	-0.0087	0.0234	0.0040	-0.0083	-0.0574	0.9563	-0.2168	-0.3662	-0.2441	2.3193	7.0101	0.0037	0.0014	-0.0040	0.7908	0.5354	-0.2427	1.4038	-0.3052	

Raw Data: S15VR23.mat
Subjects 1-5, Session VR2 & 3


29 Columns
24 time based signals
5 classification labels

Subject
Session
Segment
Arousal
Valence


344817 Rows @ 250Hz Sample Rate
23 minutes total
4m 35s per Subject
each Segment 90s

Rose and I, Rollercoaster, Swing

Data Visualization – Pendulum Swing Heart Rate



Data Visualization – Pendulum Swing EEG



Data Preparation


“Data scientists, according to interviews and expert estimates, spend from 50 percent to 80 percent of their time mired in this more mundane labor of collecting and preparing unruly digital data, before it can be explored for useful nuggets.”

For Big-Data Scientists, ‘Janitor Work’ Is Key Hurdle to Insights [New York Times August 18, 2014](#)


See also <http://blog.revolutionanalytics.com/2014/08/data-cleaning-is-a-critical-part-of-the-data-science-process.html>

Simple Segmentation and Feature Extraction

Original Data Table



Feature Table



More Complex Feature Extraction


- General-purpose signal properties
 - Time domain
 - Frequency domain
 - Other signal properties
- Caveat: these are all hand-engineered features.
- Python TSFRESH library
https://tsfresh.readthedocs.io/en/latest/text/list_of_features.html

nF.	Feature name
1	Average
2	Standard deviation
3-8	PSD- peaks frequency
9-14	PSD-peaks amplitude
15	Energy
16	Zero crossing rate
17	Energy entropy
18	Spectral centroid
19	Spectral spread
20	Spectral entropy
21	Spectral Rolloff point
22-26	MODWT- Energy of Wavelet coefficients
27-31	MODWT- Percentage of Energy of Wavelet coefficients
32-36	MODWT- Standard deviation of Wavelet coefficients
37-41	MODWT- Mean of Wavelet coefficients
42	Tsallis entropy
43	Renyi entropy
44	Shannon entropy
45-54	RSP of subbands
55	RSP- Slow wave bands-spectral bands Delta
56	RSP- Slow wave bands-spectral bands Theta
57	RSP- Slow wave bands-spectral bands Alpha
58-72	Harmonic parameters
73	Hjorth parameters– Activity
74	Hjorth parameters– Mobility
75	Hjorth parameters– Complexity
76	Skewness
77	Kurtosis
78-87	Autoregressive parameters
88-90	Percentile 25, 50, 75

Table XIV

LIST OF EXTRACTED FEATURES WITH THEIR RELEVANT NUMBERS.

The Learning Process (traditional Machine Learning)



The Learning Process (Deep Learning)







Image credit: Marco Del Pra

<https://towardsdatascience.com/time-series-classification-with-deep-learning-d238f0147d6f>

Case Study 2: EEG Biosignal Analysis

What is EEG?


- ❖ Electroencephalography (EEG) is an electrophysiological monitoring method to record electrical activities of the brain.
- ❖ It is typically noninvasive, with the electrodes placed along the scalp.
- ❖ EEG measures voltage fluctuations resulting from ionic current within the neurons of the brain.



Epileptic spike and wave discharges monitored with EEG

EEG Advantages

- ❖ The greatest advantages of EEG is its high **temporal resolution**. EEG can determine the relative strengths and positions of electrical activity in different brain regions.
- ❖ Does not require the head/body to be fixed.
- ❖ Lightweight, portable.
- ❖ Affordable.
- ❖ Easy to use.



EEG



fMRI / PET



MEG

EEG usage


- ❖ **For medical use**, EEG is one of the main diagnostic tests for epilepsy. And it might also be helpful for diagnosing or treating the following disorders:
 - Brain tumor
 - Brain damage from head injury
 - Stroke
 - Sleep disorders
 - Brain dysfunction
 - Test of brain death in comatose patients
 - Test drug effects, Etc...
- ❖ **For research use**, EEG, and the related study of ERPs are used extensively in **neuroscience, cognitive science, cognitive psychology, neurolinguistics and psychophysiological research**, etc.
- ❖ Machine learning techniques have been used on EEG data and have shown high levels of success in **classifying mental states** (Relaxed, Neutral, Concentrating), **mental emotional states** (Negative, Neutral, Positive), **Brain computer interfaces** (BCI), etc.

Challenges of EEG Processing




- Low signal-to-noise ratio (SNR). Easy to be affected by environment noise.
- Non-stationary. Statistics vary across time.
- High inter-subject variability. Physiological differences between individuals.
- A Variety of data recording settings. Different stimuli, sampling rates, number of channels.
- Require domain experts to add labels. Time consuming, data shortage.
- Domain approaches are used to pre-process data.

Is there hope for successful EEG analysis?



Deep learning (DL) could significantly simplify processing pipelines by allowing automatic end-to-end learning of preprocessing, feature extraction and classification modules, it potentially can reach competitive performance on the target task.






Deep learning as a data driven approach, it requires a large amount of data to generate considerable results comparing to traditional methods.




Many deep learning techniques have been used to improve EEG data classification accuracy but few of them tackle the issue to solve the data shortage problem in this field.


Visualizations




Common preprocessing steps



Concatenate raw data




Segment data




Amplify each channel
signal and normalize
to 0-1

Methodology

Classify with **traditional machine learning** models



Feature extractor:
e.g. Common spatial pattern (CSP)



Machine learning Method:


K Nearest Neighbors, Linear SVM, RBF SVM, Gaussian Process, Decision Tree, Random Forest, AdaBoost, Naive Bayes, Logistic Regression, LinearDiscriminantAnalysis, QuadraticDiscriminantAnalysis, etc.



Training method: `sklearn.Pipeline(CSP + MLModel)`

Methodology

Classify with deep convolutional neural networks




EEGNet

Layer (type)	Output Shape	Param #
input_2 (InputLayer)	[(None, 1, 64, 160)]	0
conv2d_2 (Conv2D)	(None, 8, 64, 160)	256
batch_normalization_1 (Batch Normalization)	(None, 8, 64, 160)	32
depthwise_conv2d (Depthwise Conv2D)	(None, 16, 1, 160)	1024
batch_normalization_2 (Batch Normalization)	(None, 16, 1, 160)	64
activation_3 (Activation)	(None, 16, 1, 160)	0
average_pooling2d_1 (Average Pooling2D)	(None, 16, 1, 40)	0
dropout_1 (Dropout)	(None, 16, 1, 40)	0
separable_conv2d (Separable Conv2D)	(None, 16, 1, 40)	512
batch_normalization_3 (Batch Normalization)	(None, 16, 1, 40)	64
activation_4 (Activation)	(None, 16, 1, 40)	0
average_pooling2d_2 (Average Pooling2D)	(None, 16, 1, 5)	0
dropout_2 (Dropout)	(None, 16, 1, 5)	0
flatten (Flatten)	(None, 80)	0
dense (Dense)	(None, 2)	162
softmax (Activation)	(None, 2)	0

Methodology

Classify with deep convolutional neural networks




ShallowConvNet

Layer (type)	Output Shape	Param #
input_1 (InputLayer)	[(None, 1, 64, 160)]	0
conv2d (Conv2D)	(None, 40, 64, 148)	560
conv2d_1 (Conv2D)	(None, 40, 1, 148)	102400
batch_normalization (BatchNo	(None, 40, 1, 148)	160
activation (Activation)	(None, 40, 1, 148)	0
average_pooling2d (AveragePo	(None, 40, 1, 17)	0
activation_1 (Activation)	(None, 40, 1, 17)	0
dropout (Dropout)	(None, 40, 1, 17)	0
flatten (Flatten)	(None, 680)	0
dense (Dense)	(None, 2)	1362
activation_2 (Activation)	(None, 2)	0

Methodology

Classify with deep convolutional neural networks



DeepConvNet

Layer (type)	Output Shape	Param #
input_3 (InputLayer)	[(None, 1, 64, 160)]	0
conv2d_3 (Conv2D)	(None, 25, 64, 156)	150
conv2d_4 (Conv2D)	(None, 25, 1, 156)	40025
batch_normalization_4 (Batch Normalization)	(None, 25, 1, 156)	100
activation_5 (Activation)	(None, 25, 1, 156)	0
max_pooling2d (MaxPooling2D)	(None, 25, 1, 78)	0
dropout_3 (Dropout)	(None, 25, 1, 78)	0
conv2d_5 (Conv2D)	(None, 50, 1, 74)	6300
batch_normalization_5 (Batch Normalization)	(None, 50, 1, 74)	200
activation_6 (Activation)	(None, 50, 1, 74)	0
max_pooling2d_1 (MaxPooling2D)	(None, 50, 1, 37)	0
dropout_4 (Dropout)	(None, 50, 1, 37)	0
conv2d_6 (Conv2D)	(None, 100, 1, 33)	25100
batch_normalization_6 (Batch Normalization)	(None, 100, 1, 33)	400
activation_7 (Activation)	(None, 100, 1, 33)	0
max_pooling2d_2 (MaxPooling2D)	(None, 100, 1, 16)	0
dropout_5 (Dropout)	(None, 100, 1, 16)	0
conv2d_7 (Conv2D)	(None, 200, 1, 12)	100200
batch_normalization_7 (Batch Normalization)	(None, 200, 1, 12)	800
activation_8 (Activation)	(None, 200, 1, 12)	0
max_pooling2d_3 (MaxPooling2D)	(None, 200, 1, 6)	0
dropout_6 (Dropout)	(None, 200, 1, 6)	0
flatten_1 (Flatten)	(None, 1200)	0
dense_1 (Dense)	(None, 2)	2402
activation_9 (Activation)	(None, 2)	0

Practical examples with source code

- Visit the shared Google Drive folder: <https://drive.google.com/drive/folders/1-GSnRyWCjodfnnkeiY5mu9pDLxAdAs2p?usp=sharing>

OR

- GitHub repository: https://github.com/imics-lab/biosignal_analysis_tutorials